



11 January 2021

Last week all X-ray, CT, Ultrasound, MRI, Nuc Med and Mammography images and reports became available in a new web viewing program. Why is this change necessary? Adobe no longer supports Flash Player (our old web product was Flash based) and Internet browsers were set to disabled Flash at the end of December.

External to Epic, the following link will still work to view images in the new program: <https://pacs.cpgph.org>

If you've previously book-marked a link for PACS please double check the link is correct. Additionally, the credentials used to log-in are what you use for CPH, not PACS and not Epic.

Getting Started

Log in to view your study exams and reports

Log In

User ID:

Password:

Preferred Language:

Enter your credentials and select a preferred language.

Using the Studies tab, enter patient search information and press Enter or click Search

Last Name: First Name: Patient ID: Date of Birth: Modality:
 Accession No.: Institution: Study Performed: Ref. Physician: Issuer of PID:
 My Patients Only

Search Preset:

Additional Search Features

Save allows you to save the current search criteria as a **Search Preset** or **My Search**.

Clear resets all fields except for the Study Performed date.

Select any **Search Preset** from the drop-down and click **Search**.

Highlight a patient study from the search results and select a study viewer to open with

or

Double-click a patient study to open it in your default study viewer

Last Name: First Name: Patient ID: Date of Birth: Modality:
 Accession No.: Institution: Study Performed: Ref. Physician: Issuer of PID:
 My Patients Only

Search Preset:

Group By	Study Date	Modality	Gender	Patient Name	Description	Patient ID	# Images	Last Name	Patient DOB	Body Part	Accession
	07-Aug-2009 8:40:13 AM	DT.DOC.MR	F	MULTIPLEXORS, JANE	MRI BREASTS BLAT	A00137	2839	MULTIPLEXORS	08-Sep-1948		A00137-3
	01-Dec-2013 12:00:00 AM	PLCT	F	MULTIPLEXORS, JANE	PET CT SKULL/THGH	A00137	882	MULTIPLEXORS	08-Sep-1948		A00137-20
	07-Apr-2012 12:00:00 AM	CT	F	MULTIPLEXORS, JANE	CTA HEAD W/ CONTRAST	A00137	775	MULTIPLEXORS	08-Sep-1948		A00137-19
	12-Jan-2018 11:10:28 AM	DOCCT	F	MULTIPLEXORS, JANE	CT RIGHT UPPER EXTREMITY W/ CONTRAST	A00137	527	MULTIPLEXORS	08-Sep-1948		A00137-9
	08-Apr-2018 4:03:15 PM	DOC.PLCT	F	MULTIPLEXORS, JANE	CT SHOULDS W/ CONTRAST	A00137	489	MULTIPLEXORS	08-Sep-1948		A00137-16

Viewer:


Click the link for "Enterprise" Viewer




From within Epic the link will take you directly to the new web viewing program without requesting log-in credentials as per current workflow. Basic functionality is as follows:

[Click here to view prior exams](#)

Enterprise Study Viewer



The yellow arrow below is pointing to the report icon - found in the upper right corner of your screen



Contact the CPH helpdesk at x4701 for assistance